Identification Data



November 16, 2020

LAB GROWN DIAMOND Certificate No: 303110197

Gemprint

Gemprint is the unique optical fingerprint for positive identification of your lab grown diamond. Register your lab grown diamond at www.Gemprint.com and receive insurance discounts up to 10%.



Laser Inscription:

The illustration depicts enlarged and approximate appearances of the inscriptions. Girdle laser inscribed "LAB GROWN" and "LG303110197"





580 Fifth Avenue LL-05, NY, NY 10036 USA • T 212.869.8985 • GCALUSA.com

ISO/IEC 17025 2017 ANAB L2177-1 Accredited Testing Lab



The 4Cs Grading Analysis

GCAL 303110197 LAB GROWN DIAMOND*

Carat Weight: 3.51

Cut: Ideal Shape: Round Brilliant Measurements: 9.67-9.73x6.01mm Hearts: Excellent Excellent Arrows: Optical Brilliance: Excellent Optical Symmetry: Excellent Polish: Excellent External Symmetry: Excellent Girdle Thickness: Medium-SI.Thick Culet Size: None

Color: E Fluorescence: None

Clarity: Identifying Characteristic(s) Characteristic Location(s): VVS1 Pinpoints Throughout Pavilion

*Comments: This man-made diamond was grown in a laboratory by the HPHT method, and has the same chemical, physical, and optical properties as a natural earth mined diamond. This diamond is Type II, which means it is devoid of nitrogen impurities.

Photomicrographs:

Actual images of the crown (top) and pavilion (bottom) of this lab grown diamond photographed at magnifications up to 10x.





Light Performance Profile

Optical Brilliance Analysis:



Brilliance is the overall return of light to the viewer. The brilliance image is a representation of (a) white areas of light return, or brilliance, and (b) dark-blue areas of light loss.









Optical Symmetry Analysis:



The colored areas of the symmetry image are indications of light handling ability, giving a visual representation of proportions and facet alignment.









Hearts and Arrows:

Precision faceting is visualized as Hearts and Arrows when brilliant cut stones are viewed in specific lighting conditions. Each pattern is the result of facet placement and alignment.





Excellent

Excellent

Proportion Diagram:

The proportion diagram illustrates the actual dimensions as recorded by optical scanning technology.

